

## Introduction


The characterization of the microscopic dynamics of functional films & coatings is an important step in every new development or quality control. We use a dynamic light scattering technique to investigate the microstructure evolution (evaporation, packing, coalescence...) of a large panoply of materials at constant or in-cresing temperature (RT-250°C) and with humidity control. This technology allows to compare different formulation and to detect the characteristic steps (Drying/curing time, phase transition...) with a very handy sampling protocol, high sensitivity and on realistic conditions (temperature, humidity, substrate...)

## The solution

- ✓ **Objectivity and accuracy to monitor and compare the curing/drying (steps, times...) of different products**
- ✓ **Possibility to analyse any sample and on any substrate**
- ✓ **Realistic conditions: Temperature & humidity control**




## How does it work [1]



## Coatings, Paints & Inks

**Characteristic times of the drying** (Open-time, Surface dry, Dry-through...)  
**Microstructure change** (packing, particles deformation, coalescence...)

Exploring temperature, humidity, substrate, thickness influence




**Clear identification of the different drying steps**

## Powder Coatings


**Clear identification of the different curing steps**

Micro-Dynamics (mD)




**Determination of the optimal curing protocol**

## Functional Film




**Clear temperature and humidity influence**  
→ **Optimize the drying protocol & process**



**Compare, rank and screen different formulations and drying conditions**

## Polymers phase transition



**Detection of polymers phase transition**

## Conclusion

Innovative **in-situ, non-invasive and handy** method to better understand your different materials allowing to:

- ✓ **Monitor and know precisely the curing and drying kinetics**
- ✓ **Determine the characteristic times of the curing process**
- ✓ **Analyse from RT up to 250°C with humidity control**
- ✓ **Evaluate the impact of the formulation, the temperature or the substrate**
- ✓ **Optimize the manufacturing protocol**

## References

- [1] D. J. Pine, D.A.W., J. X. Zhu and E. Herbolzheimer, Diffusing-wave spectroscopy: dynamic light scattering in the multiple scattering limit. Journal de physique, 1990.
- [2] Charde, S. J., Sonawane, S. S., Sonawane, S. H., & Shimpi, N. G. (2018). Degradation Kinetics of Polycarbonate Composites: Kinetic Parameters and Artificial Neural Network. Chemical and biochemical engineering quarterly, 32(2), 151-165.