

Perceptive Engineering Company

Who we are:

- 30 Employees
- Offices in the UK, Singapore, Ireland

What We Do:

- We develop software for the automation industries:
 - PAT, Advanced Process Control, Monitoring and Optimisation.

Academic and Innovation Alliances

- Universities of Cambridge, Manchester, Newcastle, Rutgers, Limerick, Strathclyde, Leeds, Surrey
- Centre for Process Innovation (CPI), Centre for Continuous Manufacturing and Crystallisation (CMAC), Institute of Chemical and Engineering Sciences (ICES Singapore), Synthesis & Solid State Pharmaceutical Centre (SSPC)
- Industrial partnerships with Siemens and GEA



Sci-Tech Daresbury UK



Singapore



The PharmaMV Platform

- Philosophy

A software platform for Advanced Process Control in the Pharmaceutical Industry: Some of the challenges....

- Management of "Instrument and Data" Integrity
- Capability to deal with Batch and Continuous Processes
- Interconnectivity Use of existing industrial standards.
- Delivering Fault tolerant Monitoring, Control, Optimization.
- Traceable User Actions and support for 21CFR Part 11.

Blending the Pharmaceutical Scientist and Automation Engineering disciplines



Perceptive Engineering

www.PerceptiveAPC.com

PharmaMV - Key Benefits Process Modelling, Monitoring and Control



		COC_Overnme_2			
Dir	rect Compr	ession Monitoring an	d Control Das	hboar	d
227.90	1 2 3	55art 10 May 10 02 51 10 May 10 02 1	52 13.130 Range 1m 5	SQ Da Dry	e 13 May 10 02 *
4.18		100 160 161 150			
let Thickness FV	And the second se		•••	•••	••••
		13 Max 10 0	9.52.00	10 10	ov 18 00 51 00 D May 18 00 51 2
ART Hardness-Pr	18 May 18 53:05:45				
Ref Hardness-Pr	18 May 18 53:08 45 Most Recent Una	knowledged Events since: 2018-05-17	09:49:14.191		
Aat Hardnass-Py	TB May 10 03:00 45 Most Recent Use	knowledged Events since: 2018-05-17 Message	09:49:14.191 Timestang	User	Acknowledged
BB 55	To May 10 03:00-45 Most Recent Une Event Id 12	knowledged Events since: 2018-05-17 Message Event Trippered 12 EV Event 12 1 Tabletikelete o mide events 1 junit	09.49-14.191 Timestamp 2018-05-18.09.52.58 2018-05-18.09.52.58	User system	Acknowledged

Advanced Process Control

- Model Predictive Control of Critical Quality Attributes using PAT derived or soft-sensor models.
- Integration with GEA Dynamic Process Control

Process Modelling

• Creation and execution of dynamic process, chemometric, mechanistic and soft-sensor models.

Data Visualisation

- Uni-variate and multi-variate visualisation tools including SPC charts, contribution/score plots, trends
- Data visualisations can be added to web-enabled dashboards and reports.

Multi-Variate Statistical Process Control

- Multivariate process monitors using Principal Component Analysis, Partial Least Squares
- MVSPC Plots with SPE & T2 statistics provided for real time monitoring.
- Implementation of Real-time release

© Perceptive Engineering

www.PerceptiveAPC.com

Process Modelling Setting the scene for "Statistical" and "Mechanistic" Models

Statistical, Empirical or Data Driven Models

- Control: Dynamic Models
- Calibration: Static Models

Created from

- Designed Plant Tests
- Historical Process Data

Mechanistic Models

- A model is one where the basic elements of the model have a direct correspondence to the underlying mechanisms in the system being modelled (1)
- Parameterised from experimental data





(1) https://grey.colorado.edu/oreilly/index.php/Mechanistic_Model

Perceptive Engineering
www.PerceptiveAPC.com



APC Case Study: TSWG

onfidential © Perceptive Engineering

www.PerceptiveAPC.com

PerceptiveAPC Platform Case Study: Twin-Screw Wet Granulator APC

Closed loop control of Critical Quality Attributes

- Twin Screw-Wet Granulation Example (GSK)
- Model Predictive Control, PAT and Chemometrics

What are the considerations when delivering robust control?

- Uncertainty
- Data Quality Monitoring
- Measurement redundancy?





fidential © Perceptive Engineering

www.PerceptiveAPC.com



Twin Screw Wet Granulator APC Process Response Testing



Twin Screw Wet Granulator APC MPC Controller – Structure



Twin Screw Wet Granulator APC Dynamic Process Model Development



Twin Screw Wet Granulator APC Simulation and Tuning



Twin Screw Wet Granulator APC Evaluation: Closed Loop Control



Twin Screw Wet Granulator APC Evaluate the Benefits





ADDOPT Scope: Digital Design to go from Molecule to Medicine



Improve / optimise for impact





IFPAC 2018

A consortium for taking technology from medium to high TRL!



13 Feb 2018

MPC development workflow Data driven v Digital Design approaches



Data Driven Workflow



Digital Design Workflow: Mechanistic model development





61

IFPAC 2018

Statistical Model Development – Response Testing

- To identify a statistical model, Pseudo Random Binary Sequence (PRBS) step testing is applied to the flowsheet model using PharmaMV.
- The screenshot below shows the step tests on the cooling rate (°C/min) and the corresponding response of C C* (C* = set point concentration).
- This data is rich enough to identify an accurate statistical model for MPC controller.





Statistical Model Development – Identification





- The statistical model is identified using the Recursive Least Squares (RLS) algorithm. The screenshot on the left shows the C-C* response to a step change in the cooling rate.
 - The crystallisation process is non-linear towards the end of the batch as there is less material to crystallise out – this is the reason why the linear model underpredicts the concentration for the last step change (black circle).



Case Study Summary

- Benefits of digital design based workflow
 - Cost-effective development of an advanced control solution.
 - Reduced experimental effort in developing process models and control strategies.
 - Reduced wastage of the API material.
 - Minimal interruption in the process production time for step-testing/model development.
 - The research in model development can ultimately be used to improve production in commercial manufacturing plants.





PROSPECT CP Twin Screw Wet Granulation Pharmaceutics Formulation Platform for Process Development and Metrological Studies



PROJECT UPDATE

Test Bed for Innovation in Solid Formulation Manufacturability

- > Twin screw granulation with associated PAT tools
- Generate process understanding to support the development of new products & processes
- Establish control strategies for commercial manufacturing processes
- Assess new and existing process analytical tools in real-life manufacturing

PROGRESS IN 2017 AND PLANS FOR 2018

- Foundational capability being built, in place Q1-Q2 2018, fully digitally-enabled
- Commissioning Q2 2018, learning project to be contracted to begin Q3 2018
 - Demonstrate utility of infrastructure
 - Implement emerging tools and engineering approaches Propose and verify control strategy for exemplar product





Linking inprocess analysis to finished product attributes

CPI – PROSPECT CP Automation Update

Automated by a fully Industry 4.0 enabled software platform for data fusion:

- Online execution of DoE, process response tests and development of empirical models.
- User-friendly workflows for model development.
- Chemometric techniques are applied to provide estimation of Critical Quality Attributes (CQAs).
- Flexible Advanced Process Control strategies to demonstrate closed loop control with feed-forward compensation for variability in raw material properties.
- Multi-Variate Analysis can be applied to detect process abnormalities to maximise process robustness.





Perceptive Engineering

www.PerceptiveAPC.com

Ongoing Research ADDoPT Case Study - TSWG



Improving Pharmaceutical processing using APC and MVA Summary

Industrial processes can be controlled by combining PAT with APC, backed up with robust data quality monitoring:

- Co-ordinated control of process actuators (CPPs)
- Minimise the variability in product quality (moisture, granule size distribution and content uniformity)
- Respond to variations in raw material and other factors
- Risk-based system design
- Back up PAT with data driven soft-sensors for CQA prediction

Digital Design approaches can be used :

- Mechanistic models can be built using (very) small-scale experiments and then scaled up.
- Reduced wastage of the API material.
- Controller robustness and process non-linearities can be explored
- Hybrid models can be used as soft-sensors
- Use a statistical model to "soak up" residual errors



Thanks for Listening

- Acknowledgements

PSE

- Sean Birmingham
- Niall Mitchell

University of Leeds

• Prof. Kevin Roberts

CPI

- Mark Taylor
- Dave Berry
- Sofia Matrali
- Tim Addison

GSK

- Don Clancy
- Benoit Inge

Perceptive Engineering:

- David Lovett
- Mihai Rascu
- Furqan Tahir
- Eduardo Lopez-Montero

ADDoPT Partners:



